# EAR THERMOMETERS with DR. TEMP





Infrared ear thermometer

- 2 seconds measurement time
- 2 person memory -30 each with temperature reading record
- 2 dual-color arrow-shaped LEDs to indicate comparison result
- Instant comparison between the current temperature and the last stored records

- Convert ear temperature to oral, rectal or axillary temperature
- Lens filter quick release function (AV07044 only)
- Automatic power-off saves on battery life
- Reading display in Celsius or Fahrenheit
- US Patent No.: 7,854,550
- With CE 0123



#### **More Convenient!**

For AV07043 & AV07044, four temperature conversion icons are used to display the relevant temperature readings of various bodily sites/locations.




## More Hygienic!

AV07044 also has a "Lens filter quick release" function that helps to remove the used lens filter from the tip by just pressing the button.



### **Keep Tracking!**

Providing 2 user IDs, with 9 memories per each user ID, Users can keep track of the temperature change and review different records of different users when necessary.



Copyright © AVITON CARE LTD. All rights reserved.



## Instant Comparison!

Dr. Temp enables instant comparison between the current temperature and the last stored records by the indication of two dual-color arrow-shaped LEDs.